



Applied Biological Materials Inc

Tel: 1-866-757-2414
Email: info@abmGood.com
Website: www.abmGood.com

Precision™ Taq DNA Polymerase

Store at -20°C

Cat. No.	Description	Concentration	Quantity
G078	Precision™ Taq DNA Polymerase	5 U/μl	500 U

Product Description

Precision™ Taq DNA Polymerase is a thermostable enzyme that replicates DNA at 70°C to 80°C. It catalyzes the polymerization of nucleotides into duplex DNA in the 5' to 3' direction in the presence of magnesium. The enzyme has a molecular weight of approximately 90 kDa. Unlike Taq DNA Polymerase, Precision™ Taq DNA Polymerase exhibits 3' to 5' exonuclease (proofreading) activity, that enables the polymerase to correct nucleotide incorporation errors. The error rate of Precision™ Taq DNA Polymerase in PCR is 2.6×10^{-6} per nucleotide per cycle. The enzyme can be used in PCR applications that demand high fidelity and can amplify templates up to 6 kb in length. The amplified PCR products are blunt-ended.

Product Components	500 U
Precision™ Taq DNA Polymerase (5 U/μl)	100 μl
5X PCR buffer, with Mg ²⁺	2 ml
25 mM MgSO ₄	1 ml
5X GC Enhancer	1 ml

Storage Buffer Components

50 mM Tris-HCl (pH 8.2 at 25°C), 0.1 mM EDTA, 1 mM DTT, 0.05 % CHAPS and 50 % glycerol.

Unit Definition

One unit of enzyme catalyzes the incorporation of 10 nmol of deoxyribonucleotides into a polynucleotide fraction in 30 mins at 70°C.

Shipping and Storage

Upon arrival, Precision™ Taq DNA Polymerase should be stored at -20°C. Avoid repeated freeze-thaw cycles of all Precision™ Taq components to retain maximum performance. All Precision™ Taq components are stable for 1 year from the date of shipping if stored and handled properly.

Protocol

The following basic protocol serves as a general guideline and starting point for any PCR amplification. Optimal reaction conditions (incubation times, temperatures, concentration of Taq DNA Polymerase, primers, MgSO₄ and template DNA) may vary and need to be optimized for each specific PCR.

All PCR experiments should be assembled in a nuclease-free environment. In addition, DNA sample preparation, reaction set-up and subsequent reaction(s) should be performed in separate areas to avoid cross contamination.

A negative control reaction (omitting template DNA) should always be performed in tandem with sample PCR to confirm the absence of DNA contamination.

1. Add the following components to a sterile 0.2 ml PCR tube sitting on ice.

Components	Volume	Final Concentration
Template DNA	~100 ng	~2 ng/μl
Forward primer (10 μM)	1 - 2.5 μl	200 - 500 nM
Reverse primer (10 μM)	1 - 2.5 μl	200 - 500 nM
5X PCR buffer, with Mg ²⁺	10 μl	1X
25 mM MgSO ₄ (Optional)*	0 - 3 μl	1.5 - 3 mM
dNTP Mix (10 mM)	1 μl	200 μM
Precision™ Taq DNA Polymerase (5 U/μl)	0.5 - 1 μl	2.5 - 5 U
5X GC Enhancer (Optional)**	10 μl	1X
Nuclease-free H ₂ O	up to 50 μl	-

* Optimal Mg²⁺ concentration is specific to each DNA template-primer set and can only be determined experimentally.

** 5X GC Enhancer is **optional** and is only recommended for PCR amplification of GC-rich DNA templates.

• We recommend preparing a mastermix for multiple reactions to minimize reagent loss and enable accurate pipetting.

- Mix contents of tube and centrifuge briefly.
- Incubate tube in a thermal cycler at 94°C for 3 mins to completely denature the template.
- Perform 30 - 40 cycles of PCR amplification as follows:

Denature: 94°C for 30 sec

Anneal: 45 - 72°C for 30 sec

Extend: 72°C for 1 min/1 kb template

- Incubate for an additional 5 mins at 72°C and maintain the reaction at 4°C. The samples can be stored at -20°C until use.
- Analyze the amplification products by agarose gel electrophoresis and visualize by ethidium bromide or SafeView™ (Cat No. G108) staining. Use appropriate molecular weight standards.

For laboratory research only. Not for clinical applications.
For technical questions, please email us at technical@abmgood.com
or visit our website at www.abmGood.com