

# SPECTRAFUGE 24D™ DIGITAL MICROCENTRIFUGE

## Spectrafuge 24D

- 24 place solid aluminum rotor
- 16,300 x g (13,300 rpm)
- Compact design
- Full digital control
- Set, display, and convert speed in RPM and/or RCF
- Choice of 5 colours
- Available StripSpin Adapter for PCR tubes and strips
- One year warranty

Labnet  
Labnet International, Inc.



Max. Speed: 13,300 rpm  
G-force: 16,300 x g  
Capacity: 24 x 1.5ml  
Timer: 1 to 30 min/cont.  
Control: Digital  
Noise Level: 55 dB  
Motor: Brushless  
Dimensions: 9.25 x 11.5 x 8.5 in  
Weight: 18 lb.

#NLC2400-\* Spectrafuge 24D with rotor \$1695.00

\*Specify Colour

(G) Grey, (B) Blue, (P) Purple, (R) Red, or (T) Teal

#NLC2400-SS StripSpin Adapter for 0.2ml tubes \$90.00



#NLC2400-SS

**\$1695.00**  
c/w 24 Place  
Rotor and Lid

Prices Valid to March 31, 2008

**DiaMed**  
LAB SUPPLIES INC.



Customer Service: 1-800-4-DIAMED

(905) 625-6021 or 1-800-434-2633

e-mail: [diamed@diamed.ca](mailto:diamed@diamed.ca) [www.diamed.ca](http://www.diamed.ca)

Prices Subject to Change without Notice

G.S.T. & P.S.T. Extra

LabNet24D-Dec2007.QXP