



5ml Macro Tubes

Snap Cap or Screw Cap

- Locking Snap Cap or Screw Cap
- Available Sterile




Made of Acetal plastic
Reusable and freezable (-90 °C)
#DLAU660AD15 & #DLAU660AD50



#DLAU660-25B

Designed to be used
between -90 °C and +80 °C.



#DLAU900-25B

Autoclavable and freezable (-90 °C).



Cap Lock Green #DLAU660-75CLG

- RNase, DNase, Pyrogen and DNA Free
- Manufactured of exceptionally high-quality, transparent polypropylene without the use of slip agents, plasticizers, mold release agents and biocide substances
- Writing surface on flat cap and side wall for quick sample identification
- Exceptional centrifugation stability of up to 25,000 x g
- Autoclavable when open (121 °C, 20 min)
- Withstands temperatures from +121 °C to -196 °C (not to be immersed in LN2 liquid phase)

Cat. #	Description	Qty.	Price
DLAU660-75N	5ml Locking Snap Cap, Graduated, Natural	200	\$59.70
DLAU660-75NS	5ml Locking Snap Cap, Graduated, Natural, Sterile	200	\$59.70
DLAU660-76N	5ml Locking Snap Cap, Graduated, Natural, 2D Enable	200	\$65.68
DLAU660-76NS	5ml Locking Snap Cap, Graduated, Natural, 2D Enable, Sterile	200	\$72.15
DLAU730-5NLSS	Screw Cap Macro Tube with Cap, Natural Sterile	200	\$74.34
DLAU660-75CLG	CapLock Clip for 5.0 ml Snap Cap Tube - Green	100	\$39.77
DLAU660AD15	Universal adapters for 15ml rotors	8	\$25.86
DLAU660AD50	Universal adapters for 50ml rotors	8	\$32.25
DLAU900-25B	Rack 25 Place for 5 ml tubes - Blue	ea.	\$8.77
DLAU900-25O	Rack 25 Place for 5 ml tubes - Orange	ea.	\$8.77
DLAU900-25R	Rack 25 Place for 5 ml tubes - Red	ea.	\$8.77
DLAU660-25B	25 Place Storage box with Clear Lid - Blue	ea.	\$10.45
DLAU660-25G	25 Place Storage box with Clear Lid - Green	ea.	\$10.45
DLAU660-25P	25 Place Storage box with Clear Lid - Pink	ea.	\$10.45
DLAU660-25Y	25 Place Storage box with Clear Lid - Yellow	ea.	\$10.45

Offers Valid to June 29, 2018