

Tubes & Closures

SAVE on
TEC Series
Screwcap
Microtubes



See Page 3 for Details

SAVE on
DiaTEC
Snap Cap
Microtubes

See Page 2 for Details



FREE with Purchase*

Fidget Spinner or Cube with any order of \$200 - \$499

Box of Truffles with any order of \$500 - \$799

Flash Drive & Pen with any order of \$800 & Over

SAVE on
5ml Macro Tubes

See Page 2 for Details



Tubes & Closures FLYER 2018 OFFERS VALID UNTIL June 29, 2018
Promo Code Q12-2018 Required at Time of Order

Promotional Giveaways will be withheld from certain companies/organizations due to their policies.
Diamed reserves the right to withhold these promotional items if the orders are not going direct to the end user.

*All promo free items available only with purchase
of consumable products listed in Flyer.

TUBE PRODUCTS

Mix'n Match

Save 10% when you Buy 10 or More DiaTEC Snap Cap Microtubes

DiaTEC

Snap Cap Microtubes

- Graduated
- Rated to 35,000 x g

Code	Volume	Seal & Cap	Qty	Price
DIATEC610-3000	0.5ml	Boilproof, Standard	1000/pk	\$39.08
DIATEC610-3001	0.5ml	Low Retention, Standard	500/pk	\$30.56
DIATEC610-3167	1.5ml	Boilproof, Standard	500/pk	\$22.51
DIATEC610-3169	1.5ml	Low Retention, Standard	2x250/pk	\$34.90
DIATEC610-3170	1.5ml	Boilproof <i>Assorted Colours</i>	500/pk	\$23.50
DIATEC610-3168	1.5ml	Large, Easy Open	500/pk	\$22.51
DIATEC610-3334	2.0ml	Boilproof, Standard	500/pk	\$25.56
DIATEC610-3335	2.0ml	Low Retention, Standard	2x250/pk	\$44.28

Accessories

DIATEC610-1551	1.5 ml Microtube Pestle, Non Sterile	10/pk	\$10.70
DIATEC610-1601	1.5 - 2.0 ml Microtube Cap Lock,	100/pk	\$44.85

Premium Series Microtubes

- Graduated
- Rated to 24,000 x g

Code	Volume	Seal & Cap	Qty	Price
PRE050-N	0.5ml	Boilproof, Standard	1000/pk	\$41.86
PRE050-A	0.5ml	As Above ASSORTED	1000/pk	\$59.80
PRE150-N	1.7ml	Boilproof, Standard	1000/pk	\$39.99
PRE150-A	1.7ml	As Above ASSORTED	1000/pk	\$55.23
PRE2000S-N	2.0ml	Low Retention, Standard	200/pk	\$22.06

Advantech Series Microtubes

- Graduated
- Rated to 20,000 x g

Code	Volume	Seal & Cap	Qty	Price
AD150-N	1.5ml	Boilproof, Standard	1000/pk	\$39.99
AD150-N500	1.5ml	Boilproof, Standard	500/pk	\$22.91
AD151-N	1.5ml	Large Cap	1000/pk	\$39.99
AD151-N500	1.5ml	Large Cap	500/pk	\$22.92
AD151-A	1.5ml	Large Cap, Assorted	1000/pk	\$47.86

5ml Macro Tubes & Accessories

- Dual Graduations – at every 0.2 & 0.5 ml
- Rated to 25,000 x g*

*When properly supported in appropriate rotor.

5ml Snap Cap - Autoclavable (121°C)

		Price
MTCC2500	Natural, Bulk, 4 x 50/pk	\$69.94
MTCC2500-AST	Rainbow Pack, Green, Blue, Red, Yellow, 4 x 50/pk	\$69.94
MTCC2520	STERILE, Natural, 10 x 20/pk	\$79.65

5ml Screw Cap - Freezer Safe (-86°C)

MTCC2530	Natural, Bulk, 5 x 100/pk	\$192.74
MTCC2540	STERILE, in Foam Racks, 10 x 50/pk	\$216.34

Racks, Workstation & Adaptors

MTCC2581	25 Place Freezer Rack, Blue	\$27.14
MTCC2592	50 Place Workstation	\$77.69
MTCC2570	Adapter, 5ml in 15ml Cavity (4/pk.)	\$76.70



STERILE Microtubes

- Natural
- Graduated
- 10 x 100 / Pk.
- Rated to 16,000 x g

Code	Volume	Seal & Cap	Price
TECBT613-NS	0.65ml	Standard Boilproof	\$25.95
TECBT620-NS	1.5ml	Standard Boilproof	\$28.95
TECBT619-NS	2.0ml	Standard Boilproof	\$29.95



Diamed's Best Microtube

TEC Series Microtubes

- Graduated
- Rated to 18,000 x g

Code	Volume	Seal & Cap	Qty	Price
TEC050-N	0.5ml	Boilproof, Standard	1000/pk	\$34.02
TEC150-N	1.5ml	Boilproof, Standard	1000/pk	\$38.95
TEC151-N*	1.5ml	Boilproof, Standard *Rated to 30,000 x g	1000/pk	\$41.11
TEC2000-N	2.0ml	Boilproof, Standard	1000/pk	\$36.92



TUBE PRODUCTS

Screw Cap Microtubes

Tube & Cap Assemblies - Natural, Bulk 500/Pk.

- For a wide range of temperatures - from -196°C to +121°C (Suitable for liquid nitrogen vapour)
- 100 µl graduations
- Ideal for long-term storage
- Tested to over 20,000 x g

Non Sterile - 500/Pk.

DIATEC610-2511	0.5 ml, Skirted	\$88.93
DIATEC610-2550	1.5 ml, Conical Graduated	\$85.45
DIATEC610-2551	1.5 ml, Skirted	\$84.45
DIATEC610-2590	2.0 ml, Conical, Graduated	\$86.02
DIATEC610-2591	2.0 ml, Skirted, Graduated	\$85.58

Sterile - 500/Pk.

DIATEC610-2700	0.5 ml, Conical	\$138.04
DIATEC610-2701	0.5 ml, Skirted	\$138.04
DIATEC610-2740	1.5 ml, Conical Graduated	\$128.44
DIATEC610-2741	1.5 ml, Skirted	\$128.44
DIATEC610-2780	2.0 ml, Conical, Graduated	\$128.44
DIATEC610-2781	2.0 ml, Skirted, Graduated	\$128.44



Tube & Cap Assemblies - Natural, Bulk

- Temperature Range - from -40°C to +121°C
100 µl graduations
- Ideal for long-term storage
- Tested to 18,000 x g

Non Sterile - 1000/Pk.

TEC053SC-N	0.5 ml, Conical	\$166.38
TEC053SSC-N	0.5 ml, Skirted	\$166.02
TEC153SC-N	1.5ml, Conical, Graduated	\$159.99
TEC153SSC-N	1.5 ml, Skirted, Graduated	\$159.99
TEC2003SC-N	2.0 ml, Conical, Graduated	\$159.99
TEC2003SSC-N	2.0 ml, Skirted, Graduated	\$159.99

Sterile - 500/Pk.

TEC053SC-NS	0.5 ml, Conical	\$124.64
TEC053SSC-NS	0.5 ml, Skirted	\$124.85
TEC153SC-NS	1.5ml, Conical, Graduated	\$124.64
TEC153SSC-NS	1.5 ml, Skirted, Graduated	\$124.89
TEC2003SC-NS	2.0 ml, Conical, Graduated	\$124.82
TEC2003SSC-NS	2.0 ml, Skirted, Graduated	\$109.91



Mix'n Match

Save 5% when you Buy 5 or More TEC Series Non Sterile Screw Cap Microtubes

Round & Conical Bottom Tubes & Caps 1000/Pk.

E/S214-2001-010	8x50 mm, 1.6ml Round Bottom, PS	\$104.28	STK8580	12x55, 3ml Round Bottom, PS	500/Bag, 1000/Pk.	\$38.23
E/S214-2301-030	8x50 mm, 1.6ml Round Bottom, PP	\$101.83	STK8582	12x55, 3ml Round Bottom, PP	500/Bag, 1000/Pk.	\$35.45
E/S214-2002-010	10x50 mm, 2.5ml, Round Bottom, PS	\$101.80	STK8550	12x75, 5 ml Conical Bottom, PS	250/Bag, 1000/Pk.	\$49.55
DIATEC1275-PST	12x75 mm, 5ml, Round Bottom, PS	\$80.42	TEC1275-PP	12x75, 5ml Round Bottom, PP	250/Bag, 1000/Pk.	\$39.80
E/S214-2023-010	12x75 mm, 4.5ml, Conical Bottom, PS	\$75.59	STK8559	12mm, Flange Cap, Natural, LDPE	1000/Pk.	\$31.13
E/S214-2367-030	12x75 mm, 5ml, Round Bottom, PP	\$72.79	STK8565A	12x75 Plug Cap, Assorted Pack, LDPE	1000/Pk.	\$30.65
E/S214-2030-010	13x100 mm, 8ml, Round Bottom, PS	\$106.33	STK8569W	Thumb Cap for Glass & Plastic Tubes, Natural, LDPE	1000/Pk.	\$39.75
E/S214-2070-010	17x100 mm, 16ml, Round Bottom, PS	\$115.96	STK8595B	Hollow Top Plug Cap, 17mm, Blue, LDPE	1000/Pk.	\$41.15
E/S214-2380-030	17x100 mm, 16ml, Round Bottom, PP	\$106.51				
E/S214-2345-03N	16x125 mm, 19ml, Round Bottom, PP	\$218.33				



STK8565A



STK8569W

DiaTEC Round Bottom Tubes & Caps

DIATEC610-1100	12x75mm, 5ml Round Bottom, PS	500/Pk.	\$22.57
DIATEC610-1101	12x75mm, 5ml Round Bottom, PP	500/Pk.	\$19.95
DIATEC610-1102	12x75mm, 5ml Round Bottom, PP, Amber	500/Pk.	\$21.33
DIATEC610-1400	17x100mm, 16ml Round Bottom, PS	500/Pk.	\$58.44
DIATEC610-1401	17x100mm, 16ml Round Bottom, Natural, PP-CO	500/Pk.	\$43.74
DIATEC610-0252	12 mm Plug Caps, White, LDPE	1000/Pk.	\$21.55
DIATEC610-0253	Universal Caps, Fit 12, 13 & 16 mm Conical & Round Bottom Tubes, Natural, LDPE	1000/Pk.	\$34.38
DIATEC610-0500	12 mm Two Position Cap, Natural, LDPE	1000/Pk.	\$25.94



Transport Vials

Sterile - with Caps, 25/Bag

E/S222-3005-070	5ml, PE, 17 x 40 mm, with PE caps, Sterile	500/Pk.	\$203.24
E/S222-3005-080	5ml, PP, 17 x 40 mm, with PE caps, Sterile	500/Pk.	\$198.40
E/S222-3007-080	7ml, PP, 17 x 54 mm, with PE caps, Sterile	500/Pk.	\$209.43
E/S222-3010-080	10ml, PP, 17 x 70 mm, with PE caps, Sterile	500/Pk.	\$238.10

Non-Sterile - Bulk, Caps Sold Separately

E/S240-3005-03K	5 ml, PP, 17 x 40 mm, bulk NO Caps	1000/Pk.	\$92.57
E/S240-3010-03K	10 ml, PP, 17 x 70 mm, bulk NO Caps	1000/Pk.	\$139.68
E/S240-3615-03K	15 ml, PP; 17 x 100 mm, bulk NO Caps	1000/Pk.	\$170.12
E/S240-3010-03K	PE Cap Only for Transport Vials, Natural	1000/Pk.	\$139.68

Non-Sterile - with Caps

E/S240-3005-060	5 ml, PP, 17 x 40 mm, with PE caps	500/Pk.	\$217.82
E/S240-3007-060	7 ml, PP, 17 x 54 mm, with PE caps	500/Pk.	\$246.95
E/S240-3007-?60	7 ml, PP, 17 x 54 mm, with PE caps Specify Colour: B-Blue, G-Green, E-Emerald, L-Lavender, R-Red, Y-Yellow	500/Pk.	\$246.95
E/S240-3010-060	10 ml, PP transport vials; 17 x 70 mm, with PE caps	1000/Pk.	\$264.10
E/S240-3615-060	15 ml, PP; 17 x 100 mm, with PE caps	1000/Pk.	\$291.32



15ml Unthreaded Centrifuge Tubes, No Caps

E/S214-2415-010	PS, Conical	125/Box, 500/Pk.	\$125.77
E/S214-2415-01K	PS, Conical	Bulk 1000/Pk.	\$199.12
E/S214-2427-030	PP, Conical	125/Box, 500/Pk.	\$121.43
E/S214-2427-03K	PP, Conical	Bulk 1000/Pk.	\$193.31



Tube Closures - 1000/Pk.

DLAU804-13?	Vacucap Tube Closures, for 13 mm Tubes, PE Specify Colour: Natural, Blue, Green, Gray, Lavender, Red, Yellow	\$41.24
DLAU804-16?	Vacucap Tube Closures, for 16 mm Tubes, PE Specify Colour: Natural, Blue, Green, Gray, Lavender, Red, Yellow	\$41.24
DLAU806-?	Universal Tube Closures, for 12-16 mm Tubes, PE Specify Colour: Natural, Black, Blue, Green, Gray, Lavender, Red, Yellow	\$35.64
DLAU808-3?	Flanged Cap, for 12mm Tubes, PE Specify Colour: Natural, Blue, Green, Red, White, Yellow	\$48.56
DLAU808-10?	Flanged Cap, for 16mm Tubes, PE Specify Colour: Natural, Blue, Green, Red, White, Yellow	\$57.53
DLAU814-12?	Pierce-It Cap, for 13mm Tubes, Thermoplastic Elastomer Specify Colour: Black, Blue, Green, Gray, Lavender, Red, White, Yellow	\$117.25
DLAU814-16?	Pierce-It Cap, for 16mm Tubes, Thermoplastic Elastomer Specify Colour: Black, Blue, Green, Gray, Lavender, Red, White, Yellow	\$117.25



Scintillation Vials

DIATEC610-4500	5ml Natural, PP Tube, PE Cap (250/Bag)	1000/Pk.	\$73.81
DIATEC610-4550	7ml Natural, PP Tube, PE Cap (250/Bag)	1000/Pk.	\$71.59
DIATEC610-4551	7ml Natural, PP Tube, PE Cap (250/Tray)	1000/Pk.	\$123.82
DLAT414	6.5ml, 16x57mm, PE Tube, PP Cap, Bulk	1000/Pk.	\$180.65
DLAT414-5	6.5ml, 16x57mm, PP Tube, PE Cap, Bulk	1000/Pk.	\$180.65
DLAT440	20ml, PP, Cap - PE, Bulk	500/Pk.	\$152.56
DLAT440-1	20ml, PP, Cap - PE, (100/Tray)	500/Pk.	\$205.28



You Must use Promotional Code Q12-2018 to Receive Flyer Special Pricing
Ensure this appears on **ALL** Purchase Orders to Receive Special Prices, there can be No Exception **OFFERS VALID UNTIL June 29, 2018**

TOLL FREE: 1.800.434.2633
FAX: 1.800.293.2987
LOCAL: 905.625.6021
FAX: 905.625.6280
WEBSITE: www.diamed.ca
EMAIL: diamed@diamed.ca

FREE Shipping
On Consumable Orders
Over \$200.⁰⁰
(Ground Only)



© Diamed Lab Supplies Inc. All rights reserved. Diamed, DiaTec, DiaSoft, DiaStretch, DiaGrip, TeclG, Tec Choice are trademarks or registered trademarks of Diamed Lab Supplies Inc. **Terms & Conditions:** Prices and promotions may change without notice. Taxes and shipping charges, if applicable, are not shown. Diamed Lab Supplies Inc. is not responsible for pricing or other errors and reserves the right to modify or cancel orders if such errors occur. **Promotional & Free Items:** Must be requested at time of order to receive. Void where prohibited by institutional policy. Qualifying amount calculated before freight and applicable taxes. **Payment Terms:** Institutional Purchase Order, Visa, MasterCard, and American Express are accepted.