

## Ready • Strain™

### Cell Straining Kits

- Includes Strainer, 50ml Tube & Cap
- Preassembled - Ready to Use
- Open Directly Under the Hood
- Guaranteed Sterile - Gamma
- Packaged in Easy Open “Tear-Strip” Bags



#### Ready Stainer™

ReadyStrain™ is the first product of its kind, providing all the necessary parts for filtration of cell suspensions and combining them into a single, preassembled, ready to use, sterile packet. These kits eliminate the contamination risks associated with opening multiple packages, and assembling separate pieces (i.e. tubes strainers and caps).

<b>MTCC5040</b>	Cell Straining Kit, c/w Strainer, 50ml Tube & Screw Cap 40 µm, Blue, Sterile, Individually Wrapped	50/Pk.	<b>\$239.62</b>
<b>MTCC5070</b>	Cell Straining Kit, c/w Strainer, 50ml Tube & Screw Cap 70 µm, White, Sterile, Individually Wrapped	50/Pk.	
<b>MTCC5100</b>	Cell Straining Kit, c/w Strainer, 50ml Tube & Screw Cap 100 µm, Yellow, Sterile, Individually Wrapped	50/Pk.	

*Each kit contains a sterile SureStrain Strainer, pre-attached to a sterile 50mL tube, with a separate sterile screw-cap ready for assembly, processing and storage after the filtration is complete. Kits are supplied in easy open “tear-strip” bags and gamma sterilized.*



Offers Valid to June 29, 2018