

DiaSHIELD

Latex & Powder Free

Nitrile Exam Gloves



Ice Blue

3 mil thick

Finger Textured

Beaded Cuff

Single Use

Ambidextrous

100 Gloves Per Box

Chemotherapy Drug Tested

Diamed.ca

1-800-434-2633

diamed@diamed.ca

www.diamed.ca

DiaSHIELD

Medical Examination Gloves, Single Use Only & Non-Sterile

CHEMOTHERAPY DRUG PERMEATION RESISTANCE

Tested per ASTM D6978-05

Chemical Tested	Minimum Breakthrough Detection Time in Minutes
Carmustine (BCNU) (3,300 ppm or 3.3 mg/mL)	NR*
Cisplatin (1,000 ppm or 1.0 mg/mL)	> 240
Cyclophosphamide (Cytosan) (20,000 ppm or 20.0 mg/mL)	> 240
Dacarbazine (DTIC) (10,000 ppm or 10.0 mg/mL)	> 240
Doxorubicin Hydrochloride (2,000 ppm or 2.0 mg/mL)	> 240
Etoposide (Toposar) (20,000 ppm or 20.0 mg/mL)	> 240
Fluorouracil (50,000 ppm or 50.0 mg/mL)	> 240
Ifosfamide (50,000 ppm or 50.0 mg/mL)	> 240
Methotrexate (25,000 ppm or 25.0 mg/mL)	> 240
Mitomycin C (500 ppm or 0.5 mg/mL)	> 240
Mitoxantrone (2,000 ppm or 2.0 mg/mL)	> 240
Paclitaxel (Taxol) (6,000 ppm or 6.0 mg/mL)	> 240
Thiotepa (10,000 ppm or 10.0 mg/mL)	NR*
Vincristine Sulfate (1,000 ppm or 1.0 mg/mL)	> 240

* Not recommended for use with Thiotepa and Carmustine

- Physical properties comply with ASTM D 6319-19
- Meets ASTM D 5151-06 freedom from holes
- Meets ASTM D 6124-06 residual powder on medical gloves
- Tested for use with chemotherapy drugs as per ASTM D6978-05
- Manufactured for: **Diamed Lab Supplies Inc.** (ISO 13485:2016)
- Made in an ISO 13485:2016/ISO 9001:2015 certified facility

AVAILABLE SIZES

DiaSHIELD-LPF	Large Ice Blue Nitrile Exam Glove, 100/box
DiaSHIELD-MPF	Medium Ice Blue Nitrile Exam Glove, 100/box
DiaSHIELD-SPF	Small Ice Blue Nitrile Exam Glove, 100/box
DiaSHIELD-XSPF	Extra Small Ice Blue Nitrile Exam Glove, 100/box

Diamed.ca

Toll-free: 1.800.434.2633 Local: 905.625.6021
Fax: 1.800.293.2987 Fax: 905.625.6280

www.diamed.ca
diamed@diamed.ca

FREE SHIPPING
on consumable orders
over \$200.00 (ground only)

