

## BEADBUG™ MICROTUBE HOMOGENIZERS



- Powerful homogenization for isolation of DNA, RNA and Proteins
- Eliminates cross contamination
- Wide variety of beads available
- Also available in 3 or 24 position models

\* Additional bead sizes/materials available. Ask for details.

BEND1036	BeadBug 6 Microtube Homogenizer	\$3313.00
BEND1032-30	Triple-Pure High Impact zirconium beads - 3.0 mm	\$231.00
BEND1032-15	Triple-Pure High Impact zirconium beads - 1.5 mm	\$231.00
BEND1032-05	Triple-Pure High Impact zirconium beads - 0.5 mm	\$231.00



### FREE

myFuge™ Mini &  
Prefilled tube kits

- 8 x 1.5/2.0ml
- 4 x PCR Strips

**SPECIAL OFFER:** Order a BEND1036 & receive a FREE BENC1008 and FREE BEND1032-\*  
\*choose from BEND1032-05, -15 or -30

Promo Code: JUL312021-BMS1



## MC-24™ High Speed Microcentrifuge

- Full Color **Touch Screen**
- mySpin™ 1 button program recall
- COMBI-Rotor™ accepts tubes and strips
- Compact, only 9 x 12 x 7.75 in.
- Variable speed control, set in rpm or rcf



**Touch Control**

Max. Speed: **13,500rpm / 16,800xg**  
Max Capacity: **24 x 1.5/2.0ml**  
**2 x PCR Strips (16 x 0.2ml)**

BENC2417	MC-24™ Microcentrifuge with 24 place COMBI-Rotor	\$2079.00
BENC2417-R5	Optional 8 x 5ml rotor	\$1355.00



BENC2417-R5




**Combi-Rotor™**  
accepts microtubes  
and PCR tubes/strips

HERMLE  
Benchmark

Ask about our **2021**  
Hermle Centrifuge **BUNDLES**  
60 Years of centrifuge experience, with a new **SPIN**

## Analytical Balance Series DX

- Backlight graphical LCD display
- Quik-Cal™ Internal Calibration
- RS232 port for printer connection
- Readability 0.0001g (0.1mg)




BENW3130

## Compact UV Transilluminator

- Four powerful 302nm UV bulbs
- UV-blocking safety cover
- Scratch resistant viewing surface
- SmartDoc™ enclosure for smart phone imaging

### FREE



SmartGlow™  
Pre-Stain

**SPECIAL OFFER:** Order a BENE3100 & receive a FREE BENE4500-PS

Promo Code: JUL312021-BMS1



BENW3101A-120	Analytical Balance, Series DX, 120g x 0.0001g	\$1995.00
BENW3130	Serial Printer for Accuris Series Dx balances	\$874.00

BENE3100	Compact UV Transilluminator, includes UV filter cover	\$756.00
----------	---	----------



## myBLOCK™ Digital Dry Baths

- Digital temperature control (set and walk away)
- Hinged lid (improves uniformity & conserves energy)
- External (optional) temperature probe (measures actual sample temperature)
- Available with Quick-Flip universal blocks (supports all common tube sizes, PCR strips and PCR plates) Patent Pending
- Blocks sold separately



\* Additional blocks available. Ask for details.

BENBSH5002	myBlock II digital dry bath, dual chamber w/out block	\$821.00
BENBSWCMB	Quick-Flip Block, 24 x 1.5ml, 32 x 0.2ml & 14 x 0.5ml	\$233.00
BENBSW15	Block, 12 x 15ml centrifuge tubes	\$174.00
BENBSW50	Block, 5 x 50ml centrifuge tubes	\$174.00

### FREE Quick-Flip™ Block




- 24 x 1.5ml
- 32 x 0.2ml
- 14 x 0.5ml

**SPECIAL OFFER:** Order a BENBSH5002 & 1-BENBSWCMB and receive a FREE BENBSWCMB.

Promo Code: JUL312021-BMS1

## INCUS-SHAKER™ SHAKING INCUBATORS

- Extremely compact, just 11" wide
- TOUCH Screen control over time, speed and temperature
- Horizontally circular 19mm orbit for aeration & mixing
- Instantly exchange flask clamps with MAGIC Clamp™ accessories (PATENTED)
- Program STEP feature for linking consecutive programs



\* Additional MAGIC Clamps available. Ask for details.

BENH1001-M	Inc-Shaker Mini with non slip rubber mat	\$4496.00
BENH1000-MR	MAGIC Clamp universal platform, (9.5 x 11.5")	\$346.00
BENH1000-MR-125	MAGIC Clamp magnetic clamp, 125mL (max. 12)	\$46.00
BENH1000-MR-250	MAGIC Clamp magnetic clamp, 250mL (max. 8)	\$53.00
BENH1000-MR-500	MAGIC Clamp magnetic clamp, 500mL (max. 5)	\$53.00
BENH1000-MR-1000	MAGIC Clamp magnetic clamp, 1000mL (max. 2)	\$82.00

### FREE hybex™ Erlenmeyer Flask

#### hybex™ Media Bottle Starter Pack

(Includes 2x100ml, 3x250ml, 3x500ml and 2x1000ml)

BENB3000-SP



**SPECIAL OFFER:** Order a BENH1001-M & receive a FREE BENB3220-500 and FREE BENB3000-SP.

Promo Code: JUL312021-BMS1

## BEADBLASTER™ 96 BALL-MILL HOMOGENIZER FOR PLATES, TUBES AND JARS

- Wet, dry or cryogrinding
- Accepts a variety of tubes, plates and jars
- Large, digital touch screen control panel
- See-through safety lid



### FREE Microplate carriers & Microtube holders



**SPECIAL OFFER:** Order a BENIPD9600 & receive a FREE BENIPD9600-96 & BENIPD9600-96A2

Promo Code: JUL312021-BMS1

BENIPD9600 BeadBlaster 96 Ball Mill Homogenizer \$16168.00

## SMARTREADER™ 96 Microplate Absorbance Readers

- Intuitive, 7 inch color touch screen control panel
- Filter based system with wavelength range from 340 to 750nm
- Operates as a stand alone system, or connects to a PC or printer
- Absorbance range: 0.0 to 4.0000 Abs
- 8 channel vertical optical path, with zero dispersion
- Integrated, multispeed plate shaker

### FREE myTemp Mini Incubator

Promo Code: JUL312021-BMS1



Internal outlet

**SPECIAL OFFER:** Order a BENMR9600 or BENMR9600-T & receive a FREE BENH2200-H



\*Includes 405, 450, 492 and 630nm filters

BENMR9600\* SmartReader 96 Microplate Absorbance Reader \$8085.00  
BENMR9600-T\* SmartReader 96 with incubation \$8605.00